

Immunofluorescent Staining Protocol

Immunofluorescent staining with antibodies against intracellular cytokines and cell surface markers provides a high resolution method to identify the nature and frequency of cells which express a particular cytokines.

Buffers

| | |
|---------------------|--|
| Cell culture media: | IMDM containing 10% (v/v) FCS |
| FACS Buffer: | 1x PBS containing 2 mM EDTA, 2 mM NaN ₃ and 2% (v/v) FCS |
| Paraformaldehyde: | 4% (w/v) in 1x PBS |

Staining Protocol

- ✓ Harvest cells from tissue and prepare a single cell suspension. (Red blood cells may be removed by lysis or density gradient.)
- ✓ Adjust cell concentration to 1×10^7 /ml.
- ✓ Transfer 100 μ l of the cell suspension in a 96-well round bottom plate, corresponding to $\sim 1 \times 10^6$ cells.
- ✓ Finally, pellet cells by centrifugation (1500 rpm) for 5 min, remove supernatant and vortex plate
- ✓ Thoroughly resuspend cells in 100 μ l of FACS Buffer with a fluorochrome-conjugated monoclonal antibody specific for a cell surface antigen for 20-30 min, 4° in the dark.
- ✓ Wash cells 2x with 100 μ l FACS Buffer per well, pellet by centrifugation (1500 rpm), and remove supernatant.
- ✓ **Optional:** Thoroughly resuspend cells in 75 μ l of 4% PFA solution and incubate for 5 min at 4°C. Cell aggregation can be avoided by vortexing prior to the addition of the PFA solution.
- ✓ Wash cells 2x with 100 μ l FACS Buffer per well, pellet by centrifugation (1500 rpm), and remove supernatant.
- ✓ Finally, resuspend in 200 μ l FACS Buffer prior to flow cytometric analysis